



Dr. med. Edith Fässler

attending physician mbF
Intensivmedizin · Dept. I

Contact

Dr. med. Edith Fässler
Switzerland

edith.faessler@kssg.ch

Units

Intensivmedizin

Profile Function

attending physician mbF

Publications (5)

Kleger G, Borgquist O, Solar M, Cole J, Jaubert P, Wendel Garcia P, Maggiorini M, Otáhal M, Smid O, Faessler E, Schrag C, Karamasis G, Keeble T, Schmidbauer S, Düring J, Leithner C, Abed-Maillard S, Navarra L, Fässler E, TTM2 Trial Investigators, Nielsen N, Friberg H, Lange T, Langeland H, Vignon P, Cassina T, Bjørkholt Olsen R, McGuigan P, Awad A, Brunetti I, Undén J, Annborn M, Lundin A, Levis A, Haenggi M, Young P, Kirkegaard H, Saxena M, Hovdenes J, Bělohávek J, Cariou A, Oddo M, Wise M, Rylander C, Ullén S, Levin H, Jakobsen J, Lilja G, Cronberg T, Pelosi P, Storm C, Taccone F, Christensen S, Grejs A, Sweet K, Bewley J, Thomas M, Hollenberg J, Chew M, Nichol A, Erlinge D, Nordberg P, Morgan M, Eastwood G, Callaway C, Joannidis M, Dankiewicz J. Hypothermia versus Normothermia after Out-of-Hospital Cardiac Arrest. *N Engl J Med* 2021; 384:2283–2294.

Westhall E, Kjaergaard J, Kuiper M, Pellis T, Stammet P, Wanscher M, Wetterslev J, Wise M, Cronberg T, TTM-trial investigators, Fässler E, Hovdenes J, Hassager C, Gasche Y, Rossetti A, van Rootselaar A, Wesenberg Kjaer T, Horn J, Ullén S, Friberg H, Nielsen N, Rosén I, Åneman A, Erlinge D, Kleger G. Standardized EEG interpretation accurately predicts prognosis after cardiac arrest. *Neurology* 2016; 86:1482–90.

Bloch A, Schrag C, Fässler E, Kleger G. Classical skin lesions resembling infective endocarditis in a patient with an infected aortic composite graft. *Am J Respir Crit Care Med* 2014; 189:e66–7.

Nielsen N, Bugge J, Hingston C, Juffermans N, Koopmans M, Køber L, Langørgen J, Lilja G, Møller J, Rundgren M, Rylander C, Smid O, Werer C, Winkel P, Friberg H, TTM Trial Investigators, Fässler E, Brunetti I, Bro-Jeppesen J, Wetterslev J, Cronberg T, Erlinge D, Gasche Y, Hassager C, Horn J, Hovdenes J, Kjaergaard J, Kuiper M, Pellis T, Stammet P, Wanscher M, Wise M, Åneman A, Al-Subaie N, Boesgaard S, Kleger G. Targeted temperature management at 33°C versus 36°C after cardiac arrest. *N Engl J Med* 2013; 369:2197–206.

Kleger G, Fässler E. Can circuit lifetime be a quality indicator in continuous renal replacement therapy in the critically ill?. *Int J Artif Organs* 2010; 33:139–46.

Projects (0)

No results found.

Kantonsspital St.Gallen

Rorschacher Strasse 95

CH-9007 St.Gallen

T: +41 71 494 11 11

support.forschung@kssg.ch